

Isotype Control Mouse IgG2b FITC

PRODUCT INFORMATION

Clone: TG2b lgG2b

Catalog#: A2242/A2252

Contents: FITC conjugated antibody buffered in 10mM PBS (pH 7.0) with 0.05%

NaN₃ and 1% BSA.

DESCRIPTION

This immunoglobulin is obtained from Balb/C mouse TG2b cell line, a hapten not expressed on human cells or human cell lines. This reagent can be used as a mouse IgG2b isotype control.

PREPARATION

The immunoglobulin is purified from ascites by protein G affinity chromatography and is conjugated with FITC under optimum conditions.

USAGE

The FITC conjugation is tested for flow cytometric analysis using 20μ I/ 10^6 cells or 100μ I peripheral blood cells.

STORAGE

Store at 4°C, should not be frozen and avoid prolonged exposure to light.